

Material Compatibility Status

Exergen		
Equipment	Metrex Product	Compatibility Status
TAT5000 Temporal Thermometer	CaviCide™	No information
TAT5000 Temporal Thermometer	CaviWipes™	No information
TAT5000 Temporal Thermometer	Envirocide™	No information
TAT5000 Temporal Thermometer	CaviCide1™	No information
TAT5000 Temporal Thermometer	CaviWipes1™	No information
TAT5000 Temporal Thermometer	CaviWipes™ 2.0	No information
TAT5000 Temporal Thermometer	CaviWipes™ HP	No information
TAT5000 Temporal Thermometer	CaviWipes™ Bleach	No information
TAT-2000C	CaviCide™	No information
TAT-2000C	CaviWipes™	No information
TAT-2000C	Envirocide™	No information
TAT-2000C	CaviCide1™	No information
TAT-2000C	CaviWipes1™	No information
TAT-2000C	CaviWipes™ 2.0	No information
TAT-2000C	CaviWipes™ HP	No information
TAT-2000C	CaviWipes™ Bleach	No information

Ultrafabrics		
Equipment	Metrex Product	Compatibility Status
Brisa	CaviCide™	Compatible
Brisa	CaviWipes™	Compatible
Brisa	Envirocide™	Compatible
Brisa	CaviCide1™	No information
Brisa	CaviWipes1™	No information
Brisa	CaviWipes™ 2.0	No information
Brisa	CaviWipes™ HP	No information
Brisa	CaviWipes™ Bleach	No information
Ultraleather	CaviCide™	Compatible
Ultraleather	CaviWipes™	Compatible
Ultraleather	Envirocide™	Compatible
Ultraleather	CaviCide1™	No information
Ultraleather	CaviCide™	Compatible
Ultraleather	CaviWipes™	Compatible
Ultraleather	Envirocide™	Compatible
Ultraleather	CaviCide1™	No information
Ultraleather	CaviWipes1™	No information
Ultraleather	CaviWipes™ 2.0	No information
Ultraleather	CaviWipes™ HP	No information
Ultraleather	CaviWipes™ Bleach	No information
Ultratech	CaviCide™	Compatible
Ultratech	CaviWipes™	Compatible
Ultratech	Envirocide™	Compatible
Ultratech	CaviCide1™	No information
Ultratech	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Ultratech	CaviWipes™ 2.0	No information
Ultratech	CaviWipes™ HP	No information
Ultratech	CaviWipes™ Bleach	No information
Ultraleather	CaviWipes1™	No information
Ultraleather	CaviWipes™ 2.0	No information
Ultraleather	CaviWipes™ HP	No information
Ultraleather	CaviWipes™ Bleach	No information
Ultratech	CaviCide™	Compatible
Ultratech	CaviWipes™	Compatible
Ultratech	Envirocide™	Compatible
Ultratech	CaviCide1™	No information
Ultratech	CaviWipes1™	No information
Ultratech	CaviWipes™ 2.0	No information
Ultratech	CaviWipes™ HP	No information
Ultratech	CaviWipes™ Bleach	No information
Brisa	CaviCide™	Compatible
Brisa	CaviWipes™	Compatible
Brisa	Envirocide™	Compatible
Brisa	CaviCide1™	No information
Brisa	CaviWipes1™	No information
Brisa	CaviWipes™ 2.0	No information
Brisa	CaviWipes™ HP	No information
Brisa	CaviWipes™ Bleach	No information