

# Material Compatibility Status

| Intuitive Surgical                            |                   |                      |
|-----------------------------------------------|-------------------|----------------------|
| Equipment                                     | Metrex Product    | Compatibility Status |
| da Vinci Acessories - Energy activation cable | CaviCide™         | Compatible           |
| da Vinci Acessories - Energy activation cable | CaviWipes™        | Compatible           |
| da Vinci Acessories - Energy activation cable | Envirocide™       | Compatible           |
| da Vinci Acessories - Energy activation cable | CaviCide1™        | No information       |
| da Vinci Acessories - Energy activation cable | CaviWipes1™       | No information       |
| da Vinci Acessories - Energy activation cable | CaviWipes™ 2.0    | No information       |
| da Vinci Acessories - Energy activation cable | CaviWipes™ HP     | No information       |
| da Vinci Acessories - Energy activation cable | CaviWipes™ Bleach | No information       |
| da Vinci Acessories - Energy activation cable | CaviCide™         | Compatible           |
| da Vinci Acessories - Energy activation cable | CaviWipes™        | Compatible           |
| da Vinci Acessories - Energy activation cable | Envirocide™       | Compatible           |
| da Vinci Acessories - Energy activation cable | CaviCide1™        | No information       |
| da Vinci Acessories - Energy activation cable | CaviWipes1™       | No information       |

| <b>Equipment</b>                              | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-----------------------------------------------|-----------------------|-----------------------------|
| da Vinci Acessories - Energy activation cable | CaviWipes™ 2.0        | No information              |
| da Vinci Acessories - Energy activation cable | CaviWipes™ HP         | No information              |
| da Vinci Acessories - Energy activation cable | CaviWipes™ Bleach     | No information              |
| daVinci S System                              | CaviCide™             | Compatible                  |
| daVinci S System                              | CaviWipes™            | Compatible                  |
| daVinci S System                              | Envirocide™           | Compatible                  |
| daVinci S System                              | CaviCide1™            | No information              |
| daVinci S System                              | CaviWipes1™           | No information              |
| daVinci S System                              | CaviWipes™ 2.0        | No information              |
| daVinci S System                              | CaviWipes™ HP         | No information              |
| daVinci S System                              | CaviWipes™ Bleach     | No information              |
| daVinci Xi System                             | CaviCide™             | Compatible                  |
| daVinci Xi System                             | CaviWipes™            | Compatible                  |
| daVinci Xi System                             | Envirocide™           | Compatible                  |
| daVinci Xi System                             | CaviCide1™            | No information              |
| daVinci Xi System                             | CaviWipes1™           | No information              |
| daVinci Xi System                             | CaviWipes™ 2.0        | No information              |
| daVinci Xi System                             | CaviWipes™ HP         | No information              |
| daVinci Xi System                             | CaviWipes™ Bleach     | No information              |