

# Material Compatibility Status

| 3M               |                   |                      |
|------------------|-------------------|----------------------|
| Equipment        | Metrex Product    | Compatibility Status |
| Model 3200       | CaviWipes1™       | No information       |
| Model 3200       | CaviWipes™ 2.0    | No information       |
| Model 3200       | CaviWipes™ HP     | No information       |
| Model 3200       | CaviWipes™ Bleach | No information       |
| Model 9681       | CaviCide™         | No information       |
| Model 9681       | CaviWipes™        | No information       |
| Model 9681       | Envirocide™       | No information       |
| Model 9681       | CaviCide1™        | No information       |
| Model 9681       | CaviWipes1™       | No information       |
| Model 9681       | CaviWipes™ 2.0    | No information       |
| Model 9681       | CaviWipes™ HP     | No information       |
| Model 9681       | CaviWipes™ Bleach | Compatible           |
| Elipar Deep cure | CaviCide™         | Compatible           |
| Elipar Deep cure | CaviWipes™        | Compatible           |
| Elipar Deep cure | Envirocide™       | Compatible           |
| Elipar Deep cure | CaviCide1™        | Compatible           |
| Elipar Deep cure | CaviWipes1™       | Compatible           |
| Elipar Deep cure | CaviWipes™ 2.0    | Compatible           |
| Elipar Deep cure | CaviWipes™ HP     | No information       |
| Elipar Deep cure | CaviWipes™ Bleach | No information       |

| <b>Equipment</b>                          | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------------------------|-----------------------|-----------------------------|
| PAPR GVP Powered Air Purifying Respirator | CaviCide™             | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviWipes™            | No information              |
| PAPR GVP Powered Air Purifying Respirator | Envirocide™           | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviCide1™            | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviWipes1™           | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviWipes™ 2.0        | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviWipes™ HP         | No information              |
| PAPR GVP Powered Air Purifying Respirator | CaviWipes™ Bleach     | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviCide™             | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviWipes™            | Compatible                  |
| Reusable Respirator FF-400 Series         | Envirocide™           | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviCide1™            | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviWipes1™           | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviWipes™ 2.0        | Compatible                  |
| Reusable Respirator FF-400 Series         | CaviWipes™ HP         | No information              |
| Reusable Respirator FF-400 Series         | CaviWipes™ Bleach     | Compatible                  |

| <b>Equipment</b>                         | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|------------------------------------------|-----------------------|-----------------------------|
| True Definition Scanner                  | CaviCide™             | Compatible                  |
| True Definition Scanner                  | CaviWipes™            | Compatible                  |
| True Definition Scanner                  | Envirocide™           | Compatible                  |
| True Definition Scanner                  | CaviCide1™            | No information              |
| M-Series Headgear                        | CaviCide™             | Not recommended for use     |
| M-Series Headgear                        | CaviWipes™            | Not recommended for use     |
| M-Series Headgear                        | Envirocide™           | Not recommended for use     |
| M-Series Headgear                        | CaviCide1™            | Not recommended for use     |
| M-Series Headgear                        | CaviWipes1™           | Not recommended for use     |
| M-Series Headgear                        | CaviWipes™ 2.0        | Not recommended for use     |
| M-Series Headgear                        | CaviWipes™ HP         | No information              |
| M-Series Headgear                        | CaviWipes™ Bleach     | Compatible                  |
| Full Facepiece Respirator<br>7000 Series | CaviCide™             | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviWipes™            | No information              |
| Full Facepiece Respirator<br>7000 Series | Envirocide™           | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviCide1™            | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviWipes1™           | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviWipes™ 2.0        | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviWipes™ HP         | No information              |
| Full Facepiece Respirator<br>7000 Series | CaviWipes™ Bleach     | Compatible                  |
| Jupiter PAPR                             | CaviCide™             | No information              |

| <b>Equipment</b>  | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------|-----------------------|-----------------------------|
| Jupiter PAPR      | CaviWipes™            | No information              |
| Jupiter PAPR      | Envirocide™           | No information              |
| Jupiter PAPR      | CaviCide1™            | No information              |
| Jupiter PAPR      | CaviWipes1™           | No information              |
| Jupiter PAPR      | CaviWipes™ 2.0        | No information              |
| Jupiter PAPR      | CaviWipes™ HP         | No information              |
| Jupiter PAPR      | CaviWipes™ Bleach     | Compatible                  |
| Master Cardiology | CaviCide™             | No information              |
| Master Cardiology | CaviWipes™            | No information              |
| Master Cardiology | Envirocide™           | No information              |
| Master Cardiology | CaviCide1™            | No information              |
| Master Cardiology | CaviWipes1™           | No information              |
| Master Cardiology | CaviWipes™ 2.0        | No information              |
| Master Cardiology | CaviWipes™ HP         | No information              |
| Master Cardiology | CaviWipes™ Bleach     | Compatible                  |
| Classic II S.E.   | CaviCide™             | No information              |
| Classic II S.E.   | CaviWipes™            | No information              |
| Classic II S.E.   | Envirocide™           | No information              |
| Classic II S.E.   | CaviCide1™            | No information              |
| Classic II S.E.   | CaviWipes1™           | No information              |
| Classic II S.E.   | CaviWipes™ 2.0        | No information              |
| Classic II S.E.   | CaviWipes™ HP         | No information              |
| Classic II S.E.   | CaviWipes™ Bleach     | Compatible                  |
| Attest 390        | CaviWipes1™           | No information              |
| Attest 390        | CaviWipes™ 2.0        | No information              |
| Attest 390        | CaviWipes™ HP         | No information              |

| <b>Equipment</b>      | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-----------------------|-----------------------|-----------------------------|
| Attest 390            | CaviWipes™ Bleach     | No information              |
| Bair Hugger Model 500 | CaviCide™             | No information              |
| Bair Hugger Model 500 | CaviWipes™            | No information              |
| Bair Hugger Model 500 | Envirocide™           | No information              |
| Bair Hugger Model 500 | CaviCide1™            | No information              |
| Bair Hugger Model 500 | CaviWipes1™           | No information              |
| Bair Hugger Model 500 | CaviWipes™ 2.0        | No information              |
| Bair Hugger Model 500 | CaviWipes™ HP         | No information              |
| Bair Hugger Model 500 | CaviWipes™ Bleach     | No information              |
| Bair Hugger Model 750 | CaviCide™             | No information              |
| Bair Hugger Model 750 | CaviWipes™            | No information              |
| Bair Hugger Model 750 | Envirocide™           | No information              |
| Bair Hugger Model 750 | CaviCide1™            | No information              |
| Bair Hugger Model 750 | CaviWipes1™           | No information              |
| Bair Hugger Model 750 | CaviWipes™ 2.0        | No information              |
| Bair Hugger Model 750 | CaviWipes™ HP         | No information              |
| Bair Hugger Model 750 | CaviWipes™ Bleach     | No information              |
| Bair Hugger Model 775 | CaviCide™             | Compatible                  |
| Bair Hugger Model 775 | CaviWipes™            | Compatible                  |
| Bair Hugger Model 775 | Envirocide™           | Compatible                  |
| Bair Hugger Model 775 | CaviCide1™            | Compatible                  |
| Bair Hugger Model 775 | CaviWipes1™           | Compatible                  |
| Bair Hugger Model 775 | CaviWipes™ 2.0        | Compatible                  |
| Bair Hugger Model 775 | CaviWipes™ HP         | Compatible                  |
| Bair Hugger Model 775 | CaviWipes™ Bleach     | Compatible                  |
| Bair Paws Model 850   | CaviCide™             | No information              |

| <b>Equipment</b>     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------|-----------------------|-----------------------------|
| Bair Paws Model 850  | CaviWipes™            | No information              |
| Bair Paws Model 850  | Envirocide™           | No information              |
| Bair Paws Model 850  | CaviCide1™            | No information              |
| Bair Paws Model 850  | CaviWipes1™           | No information              |
| Bair Paws Model 850  | CaviWipes™ 2.0        | No information              |
| Bair Paws Model 850  | CaviWipes™ HP         | No information              |
| Bair Paws Model 850  | CaviWipes™ Bleach     | Compatible                  |
| Bair Paws Model 875  | CaviCide™             | No information              |
| Bair Paws Model 875  | CaviWipes™            | No information              |
| Bair Paws Model 875  | Envirocide™           | No information              |
| Bair Paws Model 875  | CaviCide1™            | No information              |
| Classic II Pediatric | CaviWipes1™           | No information              |
| Classic II Pediatric | CaviWipes™ 2.0        | No information              |
| Classic II Pediatric | CaviWipes™ HP         | No information              |
| Classic II Pediatric | CaviWipes™ Bleach     | Compatible                  |
| Model 145 Ranger     | CaviCide™             | Compatible                  |
| Model 145 Ranger     | CaviWipes™            | Compatible                  |
| Model 145 Ranger     | Envirocide™           | Compatible                  |
| Model 145 Ranger     | CaviCide1™            | Compatible                  |
| Model 145 Ranger     | CaviWipes1™           | Compatible                  |
| Model 145 Ranger     | CaviWipes™ 2.0        | Compatible                  |
| Model 145 Ranger     | CaviWipes™ HP         | No information              |
| Model 145 Ranger     | CaviWipes™ Bleach     | Compatible                  |
| Model 245            | CaviCide™             | No information              |
| Model 245            | CaviWipes™            | No information              |
| Model 245            | Envirocide™           | No information              |

| <b>Equipment</b>        | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------|-----------------------|-----------------------------|
| Model 245               | CaviCide1™            | No information              |
| Model 245               | CaviWipes1™           | No information              |
| Model 245               | CaviWipes™ 2.0        | No information              |
| Model 245               | CaviWipes™ HP         | No information              |
| Model 245               | CaviWipes™ Bleach     | Compatible                  |
| Model 3000              | CaviCide™             | No information              |
| Model 3000              | CaviWipes™            | No information              |
| Model 3000              | Envirocide™           | No information              |
| Model 3000              | CaviCide1™            | No information              |
| Model 3000              | CaviWipes1™           | No information              |
| Model 3000              | CaviWipes™ 2.0        | No information              |
| Model 3000              | CaviWipes™ HP         | No information              |
| Model 3000              | CaviWipes™ Bleach     | No information              |
| Model 3100              | CaviCide™             | No information              |
| Model 3100              | CaviWipes™            | No information              |
| Model 3100              | Envirocide™           | No information              |
| Model 3100              | CaviCide1™            | No information              |
| Model 3100              | CaviWipes1™           | No information              |
| Model 3100              | CaviWipes™ 2.0        | No information              |
| Model 3100              | CaviWipes™ HP         | No information              |
| Model 3100              | CaviWipes™ Bleach     | No information              |
| Model 3200              | CaviCide™             | No information              |
| Model 3200              | CaviWipes™            | No information              |
| Model 3200              | Envirocide™           | No information              |
| Model 3200              | CaviCide1™            | No information              |
| True Definition Scanner | CaviWipes1™           | No information              |

| <b>Equipment</b>        | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|-------------------------|-----------------------|-----------------------------|
| True Definition Scanner | CaviWipes™ 2.0        | Compatible                  |
| True Definition Scanner | CaviWipes™ HP         | No information              |
| True Definition Scanner | CaviWipes™ Bleach     | No information              |
| Cart and Mobile Edition | CaviCide™             | Compatible                  |
| Cart and Mobile Edition | CaviWipes™            | Compatible                  |
| Cart and Mobile Edition | Envirocide™           | Compatible                  |
| Cart and Mobile Edition | CaviCide1™            | No information              |
| Cart and Mobile Edition | CaviWipes1™           | No information              |
| Cart and Mobile Edition | CaviWipes™ 2.0        | Compatible                  |
| Cart and Mobile Edition | CaviWipes™ HP         | No information              |
| Cart and Mobile Edition | CaviWipes™ Bleach     | No information              |
| Bair Paws Model 875     | CaviWipes1™           | No information              |
| Bair Paws Model 875     | CaviWipes™ 2.0        | No information              |
| Bair Paws Model 875     | CaviWipes™ HP         | No information              |
| Bair Paws Model 875     | CaviWipes™ Bleach     | Compatible                  |
| Cardiology IV           | CaviCide™             | No information              |
| Cardiology IV           | CaviWipes™            | No information              |
| Cardiology IV           | Envirocide™           | No information              |
| Cardiology IV           | CaviCide1™            | No information              |
| Cardiology IV           | CaviWipes1™           | No information              |
| Cardiology IV           | CaviWipes™ 2.0        | No information              |
| Cardiology IV           | CaviWipes™ HP         | No information              |
| Cardiology IV           | CaviWipes™ Bleach     | Compatible                  |
| Clean-Trace™ LX25       | CaviCide™             | No information              |
| Clean-Trace™ LX25       | CaviWipes™            | No information              |
| Clean-Trace™ LX25       | Envirocide™           | No information              |



| <b>Equipment</b>                | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|---------------------------------|-----------------------|-----------------------------|
| Clean-Trace™ LX25               | CaviCide1™            | No information              |
| Clean-Trace™ LX25               | CaviWipes1™           | No information              |
| Clean-Trace™ LX25               | CaviWipes™ 2.0        | No information              |
| Clean-Trace™ LX25               | CaviWipes™ HP         | No information              |
| Clean-Trace™ LX25               | CaviWipes™ Bleach     | Compatible                  |
| 6000 Series Reusable Respirator | CaviCide™             | Compatible                  |
| 6000 Series Reusable Respirator | CaviWipes™            | Compatible                  |
| 6000 Series Reusable Respirator | Envirocide™           | Compatible                  |
| 6000 Series Reusable Respirator | CaviCide1™            | Compatible                  |
| 6000 Series Reusable Respirator | CaviWipes1™           | Compatible                  |
| 6000 Series Reusable Respirator | CaviWipes™ 2.0        | Compatible                  |
| 6000 Series Reusable Respirator | CaviWipes™ HP         | No information              |
| 6000 Series Reusable Respirator | CaviWipes™ Bleach     | Compatible                  |
| Adflo                           | CaviCide™             | No information              |
| Adflo                           | CaviWipes™            | No information              |
| Adflo                           | Envirocide™           | No information              |
| Adflo                           | CaviCide1™            | No information              |
| Adflo                           | CaviWipes1™           | No information              |
| Adflo                           | CaviWipes™ 2.0        | No information              |
| Adflo                           | CaviWipes™ HP         | No information              |
| Adflo                           | CaviWipes™ Bleach     | No information              |

| <b>Equipment</b>     | <b>Metrex Product</b> | <b>Compatibility Status</b> |
|----------------------|-----------------------|-----------------------------|
| Attest 390           | CaviCide™             | No information              |
| Attest 390           | CaviWipes™            | No information              |
| Attest 390           | Envirocide™           | No information              |
| Attest 390           | CaviCide1™            | No information              |
| Classic III          | CaviCide™             | No information              |
| Classic III          | CaviWipes™            | No information              |
| Classic III          | Envirocide™           | No information              |
| Classic III          | CaviCide1™            | No information              |
| Classic III          | CaviWipes1™           | No information              |
| Classic III          | CaviWipes™ 2.0        | No information              |
| Classic III          | CaviWipes™ HP         | No information              |
| Classic III          | CaviWipes™ Bleach     | Compatible                  |
| Classic II Infant    | CaviCide™             | No information              |
| Classic II Infant    | CaviWipes™            | No information              |
| Classic II Infant    | Envirocide™           | No information              |
| Classic II Infant    | CaviCide1™            | No information              |
| Classic II Infant    | CaviWipes1™           | No information              |
| Classic II Infant    | CaviWipes™ 2.0        | No information              |
| Classic II Infant    | CaviWipes™ HP         | No information              |
| Classic II Infant    | CaviWipes™ Bleach     | Compatible                  |
| Classic II Pediatric | CaviCide™             | No information              |
| Classic II Pediatric | CaviWipes™            | No information              |
| Classic II Pediatric | Envirocide™           | No information              |
| Classic II Pediatric | CaviCide1™            | No information              |

| <b>Steris</b>    |                       |                             |
|------------------|-----------------------|-----------------------------|
| <b>Equipment</b> | <b>Metrex Product</b> | <b>Compatibility Status</b> |

|                                                    |                   |                |
|----------------------------------------------------|-------------------|----------------|
| Amsco C series                                     | CaviCide™         | No information |
| Amsco C series                                     | CaviWipes™        | No information |
| Amsco C series                                     | Envirocide™       | No information |
| Amsco C series                                     | CaviCide1™        | No information |
| Amsco C series                                     | CaviWipes1™       | No information |
| Amsco C series                                     | CaviWipes™ 2.0    | No information |
| Amsco C series                                     | CaviWipes™ HP     | No information |
| Amsco C series                                     | CaviWipes™ Bleach | No information |
| Amsco CENTURY                                      | CaviCide™         | No information |
| Amsco CENTURY                                      | CaviWipes™        | No information |
| Amsco CENTURY                                      | Envirocide™       | No information |
| Amsco CENTURY                                      | CaviCide1™        | No information |
| Amsco CENTURY                                      | CaviWipes1™       | No information |
| Amsco CENTURY                                      | CaviWipes™ 2.0    | No information |
| Amsco CENTURY                                      | CaviWipes™ HP     | No information |
| Amsco CENTURY                                      | CaviWipes™ Bleach | No information |
| Amsco Evolution HC-600                             | CaviCide™         | No information |
| Amsco Evolution HC-600                             | CaviWipes™        | No information |
| Amsco Evolution HC-600                             | Envirocide™       | No information |
| Amsco Evolution HC-600                             | CaviCide1™        | No information |
| Amsco Evolution HC-600                             | CaviWipes1™       | No information |
| Amsco Evolution HC-600                             | CaviWipes™ 2.0    | No information |
| Amsco Evolution HC-600                             | CaviWipes™ HP     | No information |
| Amsco Evolution HC-600                             | CaviWipes™ Bleach | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviCide™         | No information |

|                                                    |                   |                |
|----------------------------------------------------|-------------------|----------------|
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviWipes™        | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | Envirocide™       | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviCide1™        | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviWipes1™       | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviWipes™ 2.0    | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviWipes™ HP     | No information |
| Cart & Utensil<br>Washer/disinfector AMSCO<br>1215 | CaviWipes™ Bleach | No information |
| Eagle Series 3021                                  | CaviCide™         | No information |
| Eagle Series 3021                                  | CaviWipes™        | No information |
| Eagle Series 3021                                  | Envirocide™       | No information |
| Eagle Series 3021                                  | CaviCide1™        | No information |
| Eagle Series 3021                                  | CaviWipes1™       | No information |
| Eagle Series 3021                                  | CaviWipes™ 2.0    | No information |
| Eagle Series 3021                                  | CaviWipes™ HP     | No information |
| Eagle Series 3021                                  | CaviWipes™ Bleach | No information |
| Amsco 3085 SP                                      | CaviCide™         | No information |
| Amsco 3085 SP                                      | CaviWipes™        | No information |
| Amsco 3085 SP                                      | Envirocide™       | No information |
| Amsco 3085 SP                                      | CaviCide1™        | No information |

|                             |                   |                |
|-----------------------------|-------------------|----------------|
| Amsco 3085 SP               | CaviWipes1™       | No information |
| Amsco 3085 SP               | CaviWipes™ 2.0    | No information |
| Amsco 3085 SP               | CaviWipes™ HP     | No information |
| Amsco 3085 SP               | CaviWipes™ Bleach | No information |
| 4085 General Surgical Table | CaviCide™         | Compatible     |
| 4085 General Surgical Table | CaviWipes™        | Compatible     |
| 4085 General Surgical Table | Envirocide™       | Compatible     |
| 4085 General Surgical Table | CaviCide1™        | Compatible     |
| 4085 General Surgical Table | CaviWipes1™       | Compatible     |
| 4085 General Surgical Table | CaviWipes™ 2.0    | Compatible     |
| 4085 General Surgical Table | CaviWipes™ HP     | Compatible     |
| 4085 General Surgical Table | CaviWipes™ Bleach | No information |
| 5085 SRT                    | CaviCide™         | No information |
| 5085 SRT                    | CaviWipes™        | No information |
| 5085 SRT                    | Envirocide™       | No information |
| 5085 SRT                    | CaviCide1™        | No information |
| 5085 SRT                    | CaviWipes1™       | No information |
| 5085 SRT                    | CaviWipes™ 2.0    | No information |
| 5085 SRT                    | CaviWipes™ HP     | No information |
| 5085 SRT                    | CaviWipes™ Bleach | Compatible     |
| Amsco 110LS                 | CaviCide™         | No information |
| Amsco 110LS                 | CaviWipes™        | No information |
| Amsco 110LS                 | Envirocide™       | No information |
| Amsco 110LS                 | CaviCide1™        | No information |
| Amsco 110LS                 | CaviWipes1™       | No information |
| Amsco 110LS                 | CaviWipes™ 2.0    | No information |
| Amsco 110LS                 | CaviWipes™ HP     | No information |
| Amsco 110LS                 | CaviWipes™ Bleach | No information |

|                                                     |                   |                |
|-----------------------------------------------------|-------------------|----------------|
| Amsco 400 Series                                    | CaviCide™         | No information |
| Amsco 400 Series                                    | CaviWipes™        | No information |
| Amsco 400 Series                                    | Envirocide™       | No information |
| Amsco 400 Series                                    | CaviCide1™        | No information |
| Amsco 400 Series                                    | CaviWipes1™       | No information |
| Amsco 400 Series                                    | CaviWipes™ 2.0    | No information |
| Amsco 400 Series                                    | CaviWipes™ HP     | No information |
| Amsco 400 Series                                    | CaviWipes™ Bleach | No information |
| Harmony LA surgical Lighting & Visualization system | CaviCide™         | No information |
| Harmony LA surgical Lighting & Visualization system | CaviWipes™        | No information |
| Harmony LA surgical Lighting & Visualization system | Envirocide™       | No information |
| Harmony LA surgical Lighting & Visualization system | CaviCide1™        | No information |
| Harmony LA surgical Lighting & Visualization system | CaviWipes1™       | No information |
| Harmony LA surgical Lighting & Visualization system | CaviWipes™ 2.0    | No information |
| Harmony LA surgical Lighting & Visualization system | CaviWipes™ HP     | No information |
| Harmony LA surgical Lighting & Visualization system | CaviWipes™ Bleach | No information |
| Ot 1000 Series Orthopedic                           | CaviCide™         | No information |
| Ot 1000 Series Orthopedic                           | CaviWipes™        | No information |
| Ot 1000 Series Orthopedic                           | Envirocide™       | No information |
| Ot 1000 Series Orthopedic                           | CaviCide1™        | No information |
| Ot 1000 Series Orthopedic                           | CaviWipes1™       | No information |
| Ot 1000 Series Orthopedic                           | CaviWipes™ 2.0    | No information |
| Ot 1000 Series Orthopedic                           | CaviWipes™ HP     | No information |

|                           |                   |                |
|---------------------------|-------------------|----------------|
| Ot 1000 Series Orthopedic | CaviWipes™ Bleach | Compatible     |
| Polaris Surgical Lighting | CaviCide™         | No information |
| Polaris Surgical Lighting | CaviWipes™        | No information |
| Polaris Surgical Lighting | Envirocide™       | No information |
| Polaris Surgical Lighting | CaviCide1™        | No information |
| Polaris Surgical Lighting | CaviWipes1™       | No information |
| Polaris Surgical Lighting | CaviWipes™ 2.0    | No information |
| Polaris Surgical Lighting | CaviWipes™ HP     | No information |
| Polaris Surgical Lighting | CaviWipes™ Bleach | No information |