

Material Compatibility Status

ACIST Medical Systems		
Equipment	Metrex Product	Compatibility Status
Angiographic Contrast Injection System	CaviCide™	Compatible
Angiographic Contrast Injection System	CaviWipes™	Compatible
Angiographic Contrast Injection System	Envirocide™	Compatible
Angiographic Contrast Injection System	CaviCide1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ 2.0	Compatible
Angiographic Contrast Injection System	CaviWipes™ HP	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ Bleach	Not recommended for use