

Material Compatibility Status

ACIST Medical Systems		
Equipment	Metrex Product	Compatibility Status
Angiographic Contrast Injection System	CaviCide™	Compatible
Angiographic Contrast Injection System	CaviWipes™	Compatible
Angiographic Contrast Injection System	Envirocide™	Compatible
Angiographic Contrast Injection System	CaviCide1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes1™	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ 2.0	Compatible
Angiographic Contrast Injection System	CaviWipes™ HP	Not recommended for use
Angiographic Contrast Injection System	CaviWipes™ Bleach	Not recommended for use

Planmeca		
Equipment	Metrex Product	Compatibility Status
ProSensor sensor & cable	CaviCide™	Compatible
ProSensor sensor & cable	CaviWipes™	Compatible
ProSensor sensor & cable	Envirocide™	Compatible
ProSensor sensor & cable	CaviCide1™	No information
ProSensor sensor & cable	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
ProSensor sensor & cable	CaviWipes™ 2.0	No information
ProSensor sensor & cable	CaviWipes™ HP	No information
ProSensor sensor & cable	CaviWipes™ Bleach	No information
Chair Upholstery	CaviCide™	No information
Chair Upholstery	CaviWipes™	No information
Chair Upholstery	Envirocide™	No information
Chair Upholstery	CaviCide1™	No information
Chair Upholstery	CaviWipes1™	No information
Chair Upholstery	CaviWipes™ 2.0	No information
Chair Upholstery	CaviWipes™ HP	No information
Chair Upholstery	CaviWipes™ Bleach	No information
Emerald	CaviCide™	No information
Emerald	CaviWipes™	No information
Emerald	Envirocide™	No information
Emerald	CaviCide1™	No information
Emerald	CaviWipes1™	No information
Emerald	CaviWipes™ 2.0	No information
Emerald	CaviWipes™ HP	No information
Emerald	CaviWipes™ Bleach	No information
ProScanner 2.0 sensor & cable	CaviCide™	No information
ProScanner 2.0 sensor & cable	CaviWipes™	No information
ProScanner 2.0 sensor & cable	Envirocide™	No information
ProScanner 2.0 sensor & cable	CaviCide1™	No information
ProScanner 2.0 sensor & cable	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
ProScanner 2.0 sensor & cable	CaviWipes™ 2.0	No information
ProScanner 2.0 sensor & cable	CaviWipes™ HP	No information
ProScanner 2.0 sensor & cable	CaviWipes™ Bleach	No information
Promax	CaviCide™	No information
Promax	CaviWipes™	No information
Promax	Envirocide™	No information
Promax	CaviCide1™	No information
Promax	CaviWipes1™	No information
Promax	CaviWipes™ 2.0	No information
Promax	CaviWipes™ HP	No information
Promax	CaviWipes™ Bleach	No information
ProOne	CaviCide™	Compatible
ProOne	CaviWipes™	Compatible
ProOne	Envirocide™	Compatible
ProOne	CaviCide1™	No information
ProOne	CaviWipes1™	No information
ProOne	CaviWipes™ 2.0	No information
ProOne	CaviWipes™ HP	No information
ProOne	CaviWipes™ Bleach	No information
Plan Scan	CaviCide™	No information
Plan Scan	CaviWipes™	No information
Plan Scan	Envirocide™	No information
Plan Scan	CaviCide1™	No information
Plan Scan	CaviWipes1™	No information
Plan Scan	CaviWipes™ 2.0	No information

Equipment	Metrex Product	Compatibility Status
Plan Scan	CaviWipes™ HP	No information
Plan Scan	CaviWipes™ Bleach	No information
ProSensor HD sensor & cable	CaviCide™	No information
ProSensor HD sensor & cable	CaviWipes™	No information
ProSensor HD sensor & cable	Envirocide™	No information
ProSensor HD sensor & cable	CaviCide1™	No information
ProSensor HD sensor & cable	CaviWipes1™	No information
ProSensor HD sensor & cable	CaviWipes™ 2.0	No information
ProSensor HD sensor & cable	CaviWipes™ HP	No information
ProSensor HD sensor & cable	CaviWipes™ Bleach	No information