

# Material Compatibility Status

AgaMatrix		
Equipment	Metrex Product	Compatibility Status
Jazz	CaviCide™	Not recommended for use
Jazz	CaviWipes™	Not recommended for use
Jazz	Envirocide™	Not recommended for use
Jazz	CaviCide1™	Not recommended for use
Jazz	CaviWipes1™	Not recommended for use
Jazz	CaviWipes™ 2.0	Not recommended for use
Jazz	CaviWipes™ HP	Not recommended for use
Jazz	CaviWipes™ Bleach	Not recommended for use

BBraun		
Equipment	Metrex Product	Compatibility Status
Horizon NXT	CaviWipes1™	No information
Horizon NXT	CaviWipes™ 2.0	No information
Horizon NXT	CaviWipes™ HP	No information
Horizon NXT	CaviWipes™ Bleach	Compatible
Outlook 200	CaviCide™	No information
Outlook 200	CaviWipes™	No information
Outlook 200	Envirocide™	No information
Outlook 200	CaviCide1™	No information
Outlook 200	CaviWipes1™	No information
Outlook 200	CaviWipes™ 2.0	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Outlook 200	CaviWipes™ HP	No information
Outlook 200	CaviWipes™ Bleach	Compatible
Outlook ES-Series	CaviCide™	No information
Outlook ES-Series	CaviWipes™	No information
Outlook ES-Series	Envirocide™	No information
Outlook ES-Series	CaviCide1™	No information
Outlook ES-Series	CaviWipes1™	No information
Outlook ES-Series	CaviWipes™ 2.0	No information
Outlook ES-Series	CaviWipes™ HP	No information
Outlook ES-Series	CaviWipes™ Bleach	Compatible
Vista Pumps	CaviCide™	No information
Vista Pumps	CaviWipes™	No information
Vista Pumps	Envirocide™	No information
Vista Pumps	CaviCide1™	No information
Vista Pumps	CaviWipes1™	No information
Vista Pumps	CaviWipes™ 2.0	No information
Vista Pumps	CaviWipes™ HP	No information
Vista Pumps	CaviWipes™ Bleach	Compatible
Horizon NXT	CaviCide™	No information
Horizon NXT	CaviWipes™	No information
Horizon NXT	Envirocide™	No information
Horizon NXT	CaviCide1™	No information

### **Honeywell**

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Voyager XP 1472g	CaviCide™	No information
Voyager XP 1472g	CaviWipes™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Voyager XP 1472g	Envirocide™	No information
Voyager XP 1472g	CaviCide1™	No information
Voyager XP 1472g	CaviWipes1™	No information
Voyager XP 1472g	CaviWipes™ 2.0	No information
Voyager XP 1472g	CaviWipes™ HP	Compatible
Voyager XP 1472g	CaviWipes™ Bleach	No information
Xenon 1950g.bf	CaviCide™	No information
Xenon 1950g.bf	CaviWipes™	No information
Xenon 1950g.bf	Envirocide™	No information
Xenon 1950g.bf	CaviCide1™	No information
Xenon 1950g.bf	CaviWipes1™	No information
Xenon 1950g.bf	CaviWipes™ 2.0	No information
Xenon 1950g.bf	CaviWipes™ HP	Compatible
Xenon 1950g.bf	CaviWipes™ Bleach	No information
Xenon XP 1950g	CaviCide™	Compatible
Xenon XP 1950g	CaviWipes™	Compatible
Xenon XP 1950g	Envirocide™	Compatible
Xenon XP 1950g	CaviCide1™	No information
Xenon XP 1950g	CaviWipes1™	No information
Xenon XP 1950g	CaviWipes™ 2.0	No information
Xenon XP 1950g	CaviWipes™ HP	No information
Xenon XP 1950g	CaviWipes™ Bleach	Compatible
Voyager XP 1470g	CaviCide™	No information
Voyager XP 1470g	CaviWipes™	No information
Voyager XP 1470g	Envirocide™	No information
Voyager XP 1470g	CaviCide1™	No information
Voyager XP 1470g	CaviWipes1™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Voyager XP 1470g	CaviWipes™ 2.0	No information
Voyager XP 1470g	CaviWipes™ HP	Compatible
Voyager XP 1470g	CaviWipes™ Bleach	No information