

Material Compatibility Status

Carefusion		
Equipment	Metrex Product	Compatibility Status
Alaris EtCo2	CaviWipes1™	No information
Alaris EtCo2	CaviWipes™ 2.0	No information
Alaris EtCo2	CaviWipes™ HP	No information
Alaris EtCo2	CaviWipes™ Bleach	Compatible
Alaris Syringe Module	CaviCide™	No information
Alaris Syringe Module	CaviWipes™	No information
Alaris Syringe Module	Envirocide™	No information
Alaris Syringe Module	CaviCide1™	No information
Alaris Syringe Module	CaviWipes1™	No information
Alaris Syringe Module	CaviWipes™ 2.0	No information
Alaris Syringe Module	CaviWipes™ HP	No information
Alaris Syringe Module	CaviWipes™ Bleach	Compatible
Alaris Pump IUI Connector only	CaviCide™	No information
Alaris Pump IUI Connector only	CaviWipes™	No information
Alaris Pump IUI Connector only	Envirocide™	No information
Alaris Pump IUI Connector only	CaviCide1™	No information
Alaris Pump IUI Connector only	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
Alaris Pump IUI Connector only	CaviWipes™ 2.0	No information
Alaris Pump IUI Connector only	CaviWipes™ HP	No information
Alaris Pump IUI Connector only	CaviWipes™ Bleach	Compatible
Alaris System Signature Edition gold	CaviCide™	No information
Alaris System Signature Edition gold	CaviWipes™	No information
Alaris System Signature Edition gold	Envirocide™	No information
Alaris System Signature Edition gold	CaviCide1™	No information
Alaris System Signature Edition gold	CaviWipes1™	No information
Alaris System Signature Edition gold	CaviWipes™ 2.0	No information
Alaris System Signature Edition gold	CaviWipes™ HP	No information
Alaris System Signature Edition gold	CaviWipes™ Bleach	Compatible
Avea Ventilator Systems	CaviCide™	No information
Avea Ventilator Systems	CaviWipes™	No information
Avea Ventilator Systems	Envirocide™	No information
Avea Ventilator Systems	CaviCide1™	No information
Avea Ventilator Systems	CaviWipes1™	No information
Avea Ventilator Systems	CaviWipes™ 2.0	No information
Avea Ventilator Systems	CaviWipes™ HP	No information
Avea Ventilator Systems	CaviWipes™ Bleach	Compatible

Equipment	Metrex Product	Compatibility Status
Pyxis Supply System	CaviCide™	Compatible
Pyxis Supply System	CaviWipes™	Compatible
Pyxis Supply System	Envirocide™	Compatible
Pyxis Supply System	CaviCide1™	Compatible
Pyxis Supply System	CaviWipes1™	Compatible
Pyxis Supply System	CaviWipes™ 2.0	Compatible
Pyxis Supply System	CaviWipes™ HP	No information
Pyxis Supply System	CaviWipes™ Bleach	Compatible
Alaris PCA Module	CaviCide™	No information
Alaris PCA Module	CaviWipes™	No information
Alaris PCA Module	Envirocide™	No information
Alaris PCA Module	CaviCide1™	No information
Alaris PCA Module	CaviWipes1™	No information
Alaris PCA Module	CaviWipes™ 2.0	No information
Alaris PCA Module	CaviWipes™ HP	No information
Alaris PCA Module	CaviWipes™ Bleach	Compatible
Sp02 Module	CaviCide™	No information
Sp02 Module	CaviWipes™	No information
Sp02 Module	Envirocide™	No information
Sp02 Module	CaviCide1™	No information
Sp02 Module	CaviWipes1™	No information
Sp02 Module	CaviWipes™ 2.0	No information
Sp02 Module	CaviWipes™ HP	No information
Sp02 Module	CaviWipes™ Bleach	Compatible
Alaris EtCo2	CaviCide™	No information
Alaris EtCo2	CaviWipes™	No information

Equipment	Metrex Product	Compatibility Status
Alaris EtCo2	Envirocide™	No information
Alaris EtCo2	CaviCide1™	No information

Lumenis		
Equipment	Metrex Product	Compatibility Status
A AcuPulse Carbon Dioxide Laser	CaviCide™	No information
A AcuPulse Carbon Dioxide Laser	CaviWipes™	No information
A AcuPulse Carbon Dioxide Laser	Envirocide™	No information
A AcuPulse Carbon Dioxide Laser	CaviCide1™	No information
A AcuPulse Carbon Dioxide Laser	CaviWipes1™	No information
A AcuPulse Carbon Dioxide Laser	CaviWipes™ 2.0	No information
A AcuPulse Carbon Dioxide Laser	CaviWipes™ HP	No information
A AcuPulse Carbon Dioxide Laser	CaviWipes™ Bleach	No information

Olympus		
Equipment	Metrex Product	Compatibility Status
Olympus Flexible Endoscopes	CaviCide™	No information
Olympus Flexible Endoscopes	CaviWipes™	No information
Olympus Flexible Endoscopes	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
Olympus Flexible Endoscopes	CaviCide1™	No information
Olympus Flexible Endoscopes	CaviWipes1™	No information
Olympus Flexible Endoscopes	CaviWipes™ 2.0	No information
Olympus Flexible Endoscopes	CaviWipes™ HP	No information
Olympus Flexible Endoscopes	CaviWipes™ Bleach	No information