

Material Compatibility Status

Carestream		
Equipment	Metrex Product	Compatibility Status
CS 3600 Family	CaviCide™	Compatible
CS 3600 Family	CaviWipes™	Compatible
CS 3600 Family	Envirocide™	Compatible
CS 3600 Family	CaviCide1™	No information
CS 3700	CaviWipes1™	Compatible
CS 3700	CaviWipes™ 2.0	Compatible
CS 3700	CaviWipes™ HP	Compatible
CS 3700	CaviWipes™ Bleach	Compatible
CS 8100	CaviCide™	Compatible
CS 8100	CaviWipes™	Compatible
CS 8100	Envirocide™	Compatible
CS 8100	CaviCide1™	Compatible
CS 8100	CaviWipes1™	Compatible
CS 8100	CaviWipes™ 2.0	Compatible
CS 8100	CaviWipes™ HP	Compatible
CS 8100	CaviWipes™ Bleach	Compatible
CS 9100	CaviCide™	No information
CS 9100	CaviWipes™	No information
CS 9100	Envirocide™	No information
CS 9100	CaviCide1™	No information
CS 9100	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
CS 9100	CaviWipes™ 2.0	No information
CS 9100	CaviWipes™ HP	No information
CS 9100	CaviWipes™ Bleach	No information
CS 9600	CaviCide™	Compatible
CS 9600	CaviWipes™	Compatible
CS 9600	Envirocide™	Compatible
CS 9600	CaviCide1™	Compatible
CS 9600	CaviWipes1™	Compatible
CS 9600	CaviWipes™ 2.0	Compatible
CS 9600	CaviWipes™ HP	Compatible
CS 9600	CaviWipes™ Bleach	Compatible
RVG Series IntraOral Scanner	CaviCide™	Compatible
RVG Series IntraOral Scanner	CaviWipes™	Compatible
RVG Series IntraOral Scanner	Envirocide™	Compatible
RVG Series IntraOral Scanner	CaviCide1™	Compatible
RVG Series IntraOral Scanner	CaviWipes1™	Compatible
RVG Series IntraOral Scanner	CaviWipes™ 2.0	Compatible
RVG Series IntraOral Scanner	CaviWipes™ HP	Compatible
RVG Series IntraOral Scanner	CaviWipes™ Bleach	Compatible
DRX Equipment	CaviCide™	No information
DRX Equipment	CaviWipes™	No information
DRX Equipment	Envirocide™	No information
DRX Equipment	CaviCide1™	Compatible
DRX Equipment	CaviWipes1™	Compatible
DRX Equipment	CaviWipes™ 2.0	No information
DRX Equipment	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
DRX Equipment	CaviWipes™ Bleach	Compatible
DryView Exterior	CaviCide™	No information
DryView Exterior	CaviWipes™	No information
DryView Exterior	Envirocide™	No information
DryView Exterior	CaviCide1™	Compatible
DryView Exterior	CaviWipes1™	Compatible
DryView Exterior	CaviWipes™ 2.0	No information
DryView Exterior	CaviWipes™ HP	Compatible
DryView Exterior	CaviWipes™ Bleach	Compatible
DRX Revolution	CaviCide™	No information
DRX Revolution	CaviWipes™	No information
DRX Revolution	Envirocide™	No information
DRX Revolution	CaviCide1™	Compatible
DRX Revolution	CaviWipes1™	Compatible
DRX Revolution	CaviWipes™ 2.0	No information
DRX Revolution	CaviWipes™ HP	No information
DRX Revolution	CaviWipes™ Bleach	No information
DRX-1	CaviCide™	No information
DRX-1	CaviWipes™	No information
DRX-1	Envirocide™	No information
DRX-1	CaviCide1™	Compatible
DRX-1	CaviWipes1™	Compatible
DRX-1	CaviWipes™ 2.0	No information
DRX-1	CaviWipes™ HP	Compatible
DRX-1	CaviWipes™ Bleach	Compatible
CS 1200	CaviCide™	No information

Equipment	Metrex Product	Compatibility Status
CS 1200	CaviWipes™	No information
CS 1200	Envirocide™	No information
CS 1200	CaviCide1™	No information
CS 1200	CaviWipes1™	No information
CS 1200	CaviWipes™ 2.0	No information
CS 1200	CaviWipes™ HP	No information
CS 1200	CaviWipes™ Bleach	Compatible
CS 3600 Family	CaviWipes1™	No information
CS 3600 Family	CaviWipes™ 2.0	No information
CS 3600 Family	CaviWipes™ HP	No information
CS 3600 Family	CaviWipes™ Bleach	Compatible
CS 3800	CaviCide™	Compatible
CS 3800	CaviWipes™	Compatible
CS 3800	Envirocide™	Compatible
CS 3800	CaviCide1™	No information
CS 3800	CaviWipes1™	No information
CS 3800	CaviWipes™ 2.0	No information
CS 3800	CaviWipes™ HP	No information
CS 3800	CaviWipes™ Bleach	No information
DR	CaviCide™	No information
DR	CaviWipes™	No information
DR	Envirocide™	No information
DR	CaviCide1™	Compatible
DR	CaviWipes1™	Compatible
DR	CaviWipes™ 2.0	No information
DR	CaviWipes™ HP	Compatible

Equipment	Metrex Product	Compatibility Status
DR	CaviWipes™ Bleach	Compatible
CS 1500	CaviCide™	No information
CS 1500	CaviWipes™	No information
CS 1500	Envirocide™	No information
CS 1500	CaviCide1™	No information
CS 1500	CaviWipes1™	No information
CS 1500	CaviWipes™ 2.0	No information
CS 1500	CaviWipes™ HP	No information
CS 1500	CaviWipes™ Bleach	Compatible
CS 3600	CaviCide™	Compatible
CS 3600	CaviWipes™	Compatible
CS 3600	Envirocide™	Compatible
CS 3600	CaviCide1™	Compatible
CS 3600	CaviWipes1™	Compatible
CS 3600	CaviWipes™ 2.0	Compatible
CS 3600	CaviWipes™ HP	Compatible
CS 3600	CaviWipes™ Bleach	Compatible
CS 3700	CaviCide™	Compatible
CS 3700	CaviWipes™	Compatible
CS 3700	Envirocide™	Compatible
CS 3700	CaviCide1™	Compatible

CUBESCAN		
Equipment	Metrex Product	Compatibility Status
BioCon-500	CaviCide™	No information
BioCon-500	CaviWipes™	No information
BioCon-500	Envirocide™	No information

Equipment	Metrex Product	Compatibility Status
BioCon-500	CaviCide1™	No information
BioCon-500	CaviWipes1™	No information
BioCon-500	CaviWipes™ 2.0	No information
BioCon-500	CaviWipes™ HP	No information
BioCon-500	CaviWipes™ Bleach	No information