

Material Compatibility Status

J H Emerson		
Equipment	Metrex Product	Compatibility Status
96-SL	CaviCide™	No information
96-SL	CaviWipes™	No information
96-SL	Envirocide™	No information
96-SL	CaviCide1™	No information
96-SL	CaviWipes1™	No information
96-SL	CaviWipes™ 2.0	No information
96-SL	CaviWipes™ HP	No information
96-SL	CaviWipes™ Bleach	No information
96DL	CaviCide™	No information
96DL	CaviWipes™	No information
96DL	Envirocide™	No information
96DL	CaviCide1™	No information
96DL	CaviWipes1™	No information
96DL	CaviWipes™ 2.0	No information
96DL	CaviWipes™ HP	No information
96DL	CaviWipes™ Bleach	No information
96TDL	CaviCide™	No information
96TDL	CaviWipes™	No information
96TDL	Envirocide™	No information
96TDL	CaviCide1™	No information
96TDL	CaviWipes1™	No information

Equipment	Metrex Product	Compatibility Status
96TDL	CaviWipes™ 2.0	No information
96TDL	CaviWipes™ HP	No information
96TDL	CaviWipes™ Bleach	No information

Masimo		
Equipment	Metrex Product	Compatibility Status
ECG Sensors	CaviCide™	No information
ECG Sensors	CaviWipes™	No information
ECG Sensors	Envirocide™	No information
ECG Sensors	CaviCide1™	No information
ECG Sensors	CaviWipes1™	No information
ECG Sensors	CaviWipes™ 2.0	No information
ECG Sensors	CaviWipes™ HP	No information
ECG Sensors	CaviWipes™ Bleach	No information
Radical-7	CaviCide™	Compatible
Radical-7	CaviWipes™	Compatible
Radical-7	Envirocide™	Compatible
Radical-7	CaviCide1™	Compatible
Radical-7	CaviWipes1™	Compatible
Radical-7	CaviWipes™ 2.0	Compatible
Radical-7	CaviWipes™ HP	No information
Radical-7	CaviWipes™ Bleach	Compatible