

# Material Compatibility Status

<b>Covidien/Medtronic/Puritan Bennett</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
800 Series/840	CaviCide™	Compatible
800 Series/840	CaviWipes™	Compatible
800 Series/840	Envirocide™	Compatible
800 Series/840	CaviCide1™	No information
800 Series/840	CaviWipes1™	No information
800 Series/840	CaviWipes™ 2.0	No information
800 Series/840	CaviWipes™ HP	Compatible
800 Series/840	CaviWipes™ Bleach	Compatible

<b>Nellcor Inc</b>		
<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Portable	CaviCide™	Compatible
Portable	CaviWipes™	Compatible
Portable	Envirocide™	Compatible
Portable	CaviCide1™	Compatible
Portable	CaviWipes1™	Compatible
Portable	CaviWipes™ 2.0	Compatible
Portable	CaviWipes™ HP	No information
Portable	CaviWipes™ Bleach	Compatible

<b>Roche</b>		
--------------	--	--

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Cobas Plus	CaviWipes1™	Compatible
Cobas Plus	CaviWipes™ 2.0	Compatible
Cobas Plus	CaviWipes™ HP	No information
Cobas Plus	CaviWipes™ Bleach	Compatible
cobas u 601	CaviCide™	No information
cobas u 601	CaviWipes™	No information
cobas u 601	Envirocide™	No information
cobas u 601	CaviCide1™	No information
cobas u 601	CaviWipes1™	No information
cobas u 601	CaviWipes™ 2.0	No information
cobas u 601	CaviWipes™ HP	No information
cobas u 601	CaviWipes™ Bleach	No information
Cobas Urisys 1100	CaviCide™	No information
Cobas Urisys 1100	CaviWipes™	No information
Cobas Urisys 1100	Envirocide™	No information
Cobas Urisys 1100	CaviCide1™	No information
Cobas Urisys 1100	CaviWipes1™	No information
Cobas Urisys 1100	CaviWipes™ 2.0	No information
Cobas Urisys 1100	CaviWipes™ HP	No information
Cobas Urisys 1100	CaviWipes™ Bleach	No information
Model 9180	CaviCide™	No information
Model 9180	CaviWipes™	No information
Model 9180	Envirocide™	No information
Model 9180	CaviCide1™	No information
Model 9180	CaviWipes1™	No information
Model 9180	CaviWipes™ 2.0	No information
Model 9180	CaviWipes™ HP	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Model 9180	CaviWipes™ Bleach	No information
Multiplate Analyzer	CaviCide™	No information
Multiplate Analyzer	CaviWipes™	No information
Multiplate Analyzer	Envirocide™	No information
Multiplate Analyzer	CaviCide1™	No information
Multiplate Analyzer	CaviWipes1™	No information
Multiplate Analyzer	CaviWipes™ 2.0	No information
Multiplate Analyzer	CaviWipes™ HP	No information
Multiplate Analyzer	CaviWipes™ Bleach	No information
CoaguChek Vantus	CaviCide™	No information
CoaguChek Vantus	CaviWipes™	No information
CoaguChek Vantus	Envirocide™	No information
CoaguChek Vantus	CaviCide1™	No information
CoaguChek Vantus	CaviWipes1™	No information
CoaguChek Vantus	CaviWipes™ 2.0	No information
CoaguChek Vantus	CaviWipes™ HP	No information
CoaguChek Vantus	CaviWipes™ Bleach	No information
CoaguChek XS Plus/meter housing	CaviCide™	No information
CoaguChek XS Plus/meter housing	CaviWipes™	No information
CoaguChek XS Plus/meter housing	Envirocide™	No information
CoaguChek XS Plus/meter housing	CaviCide1™	No information
CoaguChek XS Plus/meter housing	CaviWipes1™	No information
CoaguChek XS Plus/meter housing	CaviWipes™ 2.0	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
CoaguChek XS Plus/meter housing	CaviWipes™ HP	No information
CoaguChek XS Plus/meter housing	CaviWipes™ Bleach	Compatible
CoaguChek XS Pro	CaviCide™	No information
CoaguChek XS Pro	CaviWipes™	No information
CoaguChek XS Pro	Envirocide™	No information
CoaguChek XS Pro	CaviCide1™	No information
CoaguChek XS Pro	CaviWipes1™	No information
CoaguChek XS Pro	CaviWipes™ 2.0	No information
CoaguChek XS Pro	CaviWipes™ HP	No information
CoaguChek XS Pro	CaviWipes™ Bleach	No information
Coagu-Chek XS System	CaviCide™	No information
Coagu-Chek XS System	CaviWipes™	No information
Coagu-Chek XS System	Envirocide™	No information
Coagu-Chek XS System	CaviCide1™	No information
Coagu-Chek XS System	CaviWipes1™	No information
Coagu-Chek XS System	CaviWipes™ 2.0	No information
Coagu-Chek XS System	CaviWipes™ HP	No information
Coagu-Chek XS System	CaviWipes™ Bleach	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviCide™	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes™	No information
Coagu-Chek XS System-TEST STRIP GUIDE	Envirocide™	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviCide1™	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes1™	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes™ 2.0	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes™ HP	No information
Coagu-Chek XS System-TEST STRIP GUIDE	CaviWipes™ Bleach	Compatible
Accutrend Plus	CaviWipes1™	No information
Accutrend Plus	CaviWipes™ 2.0	No information
Accutrend Plus	CaviWipes™ HP	No information
Accutrend Plus	CaviWipes™ Bleach	No information
CoaguChek InRange	CaviCide™	No information
CoaguChek InRange	CaviWipes™	No information
CoaguChek InRange	Envirocide™	No information
CoaguChek InRange	CaviCide1™	No information
CoaguChek InRange	CaviWipes1™	No information
CoaguChek InRange	CaviWipes™ 2.0	No information
CoaguChek InRange	CaviWipes™ HP	No information
CoaguChek InRange	CaviWipes™ Bleach	No information
CoaguChek Pro II	CaviCide™	Compatible
CoaguChek Pro II	CaviWipes™	Compatible
CoaguChek Pro II	Envirocide™	Compatible
CoaguChek Pro II	CaviCide1™	Compatible
CoaguChek Pro II	CaviWipes1™	Compatible
CoaguChek Pro II	CaviWipes™ 2.0	Compatible
CoaguChek Pro II	CaviWipes™ HP	No information
CoaguChek Pro II	CaviWipes™ Bleach	Compatible

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Cobas 4800 System	CaviCide™	No information
Cobas 4800 System	CaviWipes™	No information
Cobas 4800 System	Envirocide™	No information
Cobas 4800 System	CaviCide1™	No information
Cobas 4800 System	CaviWipes1™	No information
Cobas 4800 System	CaviWipes™ 2.0	No information
Cobas 4800 System	CaviWipes™ HP	No information
Cobas 4800 System	CaviWipes™ Bleach	No information
Cobas Miras systems	CaviCide™	No information
Cobas Miras systems	CaviWipes™	No information
Cobas Miras systems	Envirocide™	No information
Cobas Miras systems	CaviCide1™	No information
Cobas Miras systems	CaviWipes1™	No information
Cobas Miras systems	CaviWipes™ 2.0	No information
Cobas Miras systems	CaviWipes™ HP	No information
Cobas Miras systems	CaviWipes™ Bleach	Compatible
Cobas Plus	CaviCide™	Compatible
Cobas Plus	CaviWipes™	Compatible
Cobas Plus	Envirocide™	Compatible
Cobas Plus	CaviCide1™	Compatible
Accu-Check Inform II	CaviCide™	No information
Accu-Check Inform II	CaviWipes™	No information
Accu-Check Inform II	Envirocide™	No information
Accu-Check Inform II	CaviCide1™	No information
Accu-Check Inform II	CaviWipes1™	No information
Accu-Check Inform II	CaviWipes™ 2.0	No information
Accu-Check Inform II	CaviWipes™ HP	No information

<b>Equipment</b>	<b>Metrex Product</b>	<b>Compatibility Status</b>
Accu-Check Inform II	CaviWipes™ Bleach	No information
Accu-check Inform II Base Unit	CaviCide™	No information
Accu-check Inform II Base Unit	CaviWipes™	No information
Accu-check Inform II Base Unit	Envirocide™	No information
Accu-check Inform II Base Unit	CaviCide1™	No information
Accu-check Inform II Base Unit	CaviWipes1™	No information
Accu-check Inform II Base Unit	CaviWipes™ 2.0	No information
Accu-check Inform II Base Unit	CaviWipes™ HP	No information
Accu-check Inform II Base Unit	CaviWipes™ Bleach	No information
Accutrend Plus	CaviCide™	No information
Accutrend Plus	CaviWipes™	No information
Accutrend Plus	Envirocide™	No information
Accutrend Plus	CaviCide1™	No information