

Material Compatibility Status

SilverFox		
Equipment	Metrex Product	Compatibility Status
Massage Tables	CaviCide™	No information
Massage Tables	CaviWipes™	No information
Massage Tables	Envirocide™	No information
Massage Tables	CaviCide1™	No information
Massage Tables	CaviWipes1™	No information
Massage Tables	CaviWipes™ 2.0	No information
Massage Tables	CaviWipes™ HP	No information
Massage Tables	CaviWipes™ Bleach	No information
Therapeutic Tables	CaviCide™	No information
Therapeutic Tables	CaviWipes™	No information
Therapeutic Tables	Envirocide™	No information
Therapeutic Tables	CaviCide1™	No information
Therapeutic Tables	CaviWipes1™	No information
Therapeutic Tables	CaviWipes™ 2.0	No information
Therapeutic Tables	CaviWipes™ HP	No information
Therapeutic Tables	CaviWipes™ Bleach	No information

Technicuff		
Equipment	Metrex Product	Compatibility Status
Cuff Model UA-1745L	CaviCide™	Compatible
Cuff Model UA-1745L	CaviWipes™	Compatible

Equipment	Metrex Product	Compatibility Status
Cuff Model UA-1745L	Envirocide™	Compatible
Cuff Model UA-1745L	CaviCide1™	Compatible
Cuff Model UA-1745L	CaviWipes1™	Compatible
Cuff Model UA-1745L	CaviWipes™ 2.0	Compatible
Cuff Model UA-1745L	CaviWipes™ HP	Compatible
Cuff Model UA-1745L	CaviWipes™ Bleach	Compatible